



Phone (800) 890 5113 ext. 215 • Fax (601) 353 0411 • saferlock@att.net • www.saferwheelchairs.com

SAFER LOCK CLEANING PROCEDURES:

The material used in the Safer Automatic Wheelchair Wheel Locks is developed by one of America's top corporations, DuPont. The plastic materials used in the locks include Zytel nylon, a component that is strong, durable and slides easily against itself and metal without lubrication. Even when exposed to water, Zytel will not corrode.

All hardware is stainless steel and zinc, allowing Safer Locks to be washed using any cleaning method without the risk of rust or corrosion.