

PREVENTING WHEELCHAIR FALLS IS EASIER THAN YOU THINK

Most wheelchair falls occur simply because the manual wheel locks are not used while transferring to and from the wheelchair, causing the chair to roll away while sitting or standing. Safer Locks Automatic operation is designed specifically to prevent these types of falls from ever occurring.



BOLT-ON DESIGN

- Simple installation on many leading wheelchairs
- No special tools required
- Easily removed if needed



NORMAL OPERATION

- Wheelchair folds as normal
- User can propel wheelchair both directions when safely seated



HUB LOCKING

- No easily damaged exposed parts
- No premature wear on tires
- When engaged, wheels lock in both directions



OVERRIDE HANDLES

Allows caregiver to disengage the locks for easy transport and when rolling the wheelchair while unoccupied



MAINTENANCE

No additional maintenance or adjusting once properly installed



WARRANTY

Unmatched TWO YEAR WARRANTY from date of purchase



MADE IN THE USA

Manufactured, assembled and serviced at our location in Jackson, Mississippi

THE PROACTIVE APPROACH TO WHEELCHAIR SAFETY



Safer Locks
Automatic
Wheelchair
Wheel Locks

For additional details visit
SaferWheelchairs.com or contact
your local Safer Locks dealer below:
